




禾瑞印象診所

## 糞便採集說明書

### 免疫法糞便檢體採集步驟

#### How To Collect Stool Sample



#### 注意事項 (Precautions) :

- 按照上面圖解採集糞便檢體，完成放入夾鏈袋中並置於陰涼處。  
Collect the stool sample as shown above, put it into the zipper bag and keep in a cool place after completed.
  - 勿將管中的液體倒出、或任意於管中加水。  
Do not pour out the liquid or add water to the sampling bottle.
  - 遇痔瘡出血或女性月經期間，暫緩採便。  
Do not collect sample during bleeding hemorrhoids or menstrual period.
- ※ 請珍惜醫療資源，將採集完畢之容器儘速交回指定醫療院所或衛生所。  
Return the collected container to the medical institutions as soon as possible.

- 1、採檢管上標籤請先寫上姓名。
- 2、在馬桶近端斜面鋪上兩張衛生紙。
- 3、如廁坐姿稍往前移，糞便落在衛生紙上，就不會滾入馬桶內沾到水。
- 4、如果對 2-3 步驟沒有把握，可先在地上鋪衛生紙。
- 5、蹲在紙上，使糞便落在衛生紙上。
- 6、旋開採檢管一端綠色蓋子，取出螺旋採集棒。
- 7、使用螺旋採集棒在糞便的不同地方來回輕劃檢體，使少量的糞便能覆蓋溝槽即可。
- 8、棒子插回採檢管，將蓋子用力押緊（有「喀」一聲才有蓋緊），插入後請不要再拔開。）
- 9、採檢管放入塑膠袋內。